

Pediatrics

Specialty
Steel

Furnishings

NOVUM

Medical Products

Single Post Mayo Stand

Model MS-SP-035



Features:

- Single post stand with "U" base design provides support and versatility
- Removable stainless steel tray measures 19" x 16.5"
- Tray height adjusts from 32"- 49" with lock
- Two, 3" casters provide additional mobility
- Base Dimensions: 22" (W) x 18.5" (D)

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com