


NOVUM
Medical Products

Behavioral Health Bed





Model # NV-BHF-B1000



Features:

- Removable top for installation
- Inside angle brackets for floor mount
- Made for the comfort and security of the patients who need it most
- Help eliminate pinch points, reduce access to assembly points and has no hooks
- Use steel cam and dowels for assembly, steel-on-steel construction for the best security available
- Cam and dowel assembly allows parts to be replaced if damaged or abused
- Rounded edges and sloped tops to reduce the risk of ligature and garroting
- Ships fully assembled

Color Options



Dimensions:

37.5"W x 82.25"L x 22"H
(Fits 36"x80" twin)

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com