



Disposable Isolation Gown



Features:

- SMS 35 gsm Coated
- Level 2 Gown
- Blue (available in Yellow and White)
- One Size Fits Most
Length - 48.5"
Chest - 29.5"
Sweep - 60.5"
Sleeve - 22"
- Three (3) Layer Fabric (available in 35gsm/45gsm/60gsm) Spunbound Polypropylene, Melt Blown and Spunbound Polypropylene Non-Reinforced Chest and Sleeve
- Passed: Level 2
- Ultrasonic Welded Seams Neck Velcro Closure or tie, waist tie and elastic or knitted cuff (Variations available) Individual Polybag & 80 pieces per cart
- 100 gowns per carton