

Trendelenburg Sleeper Recliner

Model IRC-T-LO52

Features:

- Lever operated footrest
- Fingerpull operated back
- Back locks in any position from upright to a flat sleeping position
- Dual foot pedal trendelenburg release
- Seamless back
- Easy Clean out
- Contoured back with headrest and lumbar support
- Wide shoulders for comfortable patient positioning
- 3" swivel locking transport casters

Specifications:

- Dimensions: 44"H x 28"W x 34"D
- Seat: 19"H x 24"W x 22"D
- Fully Reclined Depth: 71"
- Weight Capacity: 500 lbs.
- Warranty: Limited 7 Year



Trendelenburg Sleeper Recliner

Model IRC-T-LO52



Options:


Arm Caps



None



Wood
IRC-A-WA-A



Polyester
IRC-A-PA-A

Accessories



Push Bar
IRC-A-PB-A



Self Storing
Swivel Tablet
IRC-A-ST-A



Folding
Tablet
IRC-A-FT-A



Drainage
Bag Holder
IRC-A-DB-A



Restraint
Strap Loop
IRC-A-RS-A



IV Pole and
Holder
IRC-A-IV-A

*Rolling Foot
Tray
*Moisture
Barrier
* CAL-133

IRC-A-RF-A
IRC-A-MB-A
IRC-A-TB-A