

# Fluid Warmer

EC350



## Features:

- Programmable control. Silent operation with no fans or moving parts. No lint traps required.
- Single pane, energy-efficient, e-coated glass window in door.
- WarmSafe Zone Heating: Heat is evenly balanced throughout the chamber to provide safe fluid warming.
- Durable: The affordable, rugged stainless steel cabinets are made to last and work just as hard as you do.
- Convenient: Keep in patient or procedure rooms with this countertop model.
- Capacity: 24 1-liter bottles or 24 1-liter bags
- Size: 28 in (71.1 cm) H; 18.5 in (47 cm) W; 23.5 in (59.7 cm) D
- Mix and Match: Any size for any space. Pair with a Titan Series blanket warming cabinet to complete your patient warming system. Stacks with EC250 and EC350 blanket warmers
- Organization: Standard glide out basket and adjustable shelf. Can be used on countertop or freestanding on floor.
- Adjustable temperature range of 32° - 71° C (90° - 160° F).
- Four digit L.E.D. display.

## Specifications

- Dimensions: 28”H x 18.5”W x 23.5”D
- Capacity: 24 1-liter bottles/bags
- Safety: In the event of a power failure the cabinet will remember its programming and begins to operate as before when power is restored.
- Safety: A warming shut-off system, separate from the electronic control, prevents overheating.