

## Glider Rocker 965

### Model 965



## Features:

- Auto-locking gliding mechanism
- Extended armrests with ergonomic grip
- Slip on backrest for safety
- Lateral anti-tip supports
- Rear casters to aid movement
- Multi-position locking
- Antimicrobial wood finishes and cushion covers
- Maple hardwood base
- Treated wood to inhibit the growth of bacteria, mildew, fungi, and odors
- Finished with stain, sealer, and lacquer
- Ultra-fluid and smooth gliding motion
- Lightweight and easy to move about
- Strong and durable hardwood construction
- High-density foam for greater comfort
- Padded armrests
- Fast production lead time

## Dimensions:

Seat weight capacity:	250 lb	113.40 kg
Chair - Weight:	41 lb	18.60 Kg
Width:	25 ¾"	65.40 cm
Depth:	33 ½"	85.09 cm
Height:	43"	109.22 cm
Floor-to-seat height:	18 ½"	47 cm
Seat depth:	19"	48.26 cm
Seat width:	18 ¼"	46.36 cm
Backrest height:	27"	68.58 cm