

Pediatrics

Specialty  
Steel

Furnishings

# NOVUM

Medical Products

## 4-Outlet Medical-Grade Power Strip

*Model IVSTAND-PS4*

### Features:

- Ideal for use in operating rooms, patient-care vicinities, therapy rooms, imaging environments and administrative areas
- For use on IV poles, mobile carts, crash carts and mobile workstations, when properly mounted
- 120V, 50/60 Hz power strip complies with UL 1363A specifications
- Silver ionic powder coating is 99.9% effective in inhibiting the growth of bacteria
- Patented coating helps to reduce the risk of healthcare-acquired infections in hospitals and other medical environments
- Coating meets the requirements of JIS Z 2801:2000, an international standard for evaluating efficacy in antimicrobial products
- 4 NEMA 5-15R-HG hospital-grade outlets provide consistent power
- 15-ft. power cord provides ample reach to the power source
- Rugged, all-metal housing for long life in demanding environments
- Integrated locking outlet safety covers block access to unused outlets
- Dual, resettable built-in 15A circuit breakers protect against dangerous overloads
- Switchless design prevents accidental shut-off
- Built-in mounting tabs and flanges makes easy to install on a variety of surfaces



Phone | 716.759.7200  
Toll Free | 800.274.2742  
Fax | 716.759.0700  
Email | [info@novummed.com](mailto:info@novummed.com)  
[www.novummed.com](http://www.novummed.com)



## 4-Outlet Medical-Grade Power Strip

*Model IVSTAND-PS4*



### Specifications:

- Nominal Input Voltage(s) Supported: 120V AC
- Recommended Electrical Service: 15A, 120V
- Input Plug Features: Hospital Grade Plug
- Input Cord Gauge Type: 14 gauge; SJT
- Voltage Compatibility (VAC): 120
- Input Plug Type: NEMA 5-15P-HG
- Input Cord Length (ft.): 15
- Frequency Compatibility: 50 / 60 Hz
- Output (Watts): 1800
- Output Receptacles: (4) 5-15R-HG
- Circuit Breaker (amps): 15
- Overload Protection: 2 x 15amp breaker (1 line, 1 neutral)
- Transformer Accommodation: (2) Yes
- Outlet Orientation: Single row of 4
- Form Factors Supported: IV poles, mobile carts, crash carts and mobile medical workstations
- Included Mounting Accessories: Mounting tabs and flanges
- Mounting Holes or Slots: Yes (Includes flanges)
- Slot or Hole Mount Measurement: Distance between keyhole tabs = 127mm; end tabs = 207.45mm
- Material of Construction: Metal
- Outlets Measurement: 32.5mm
- Unit Dimensions (hwd / in.): 8.65 x 2.5 x 1.5
- Unit Weight (lbs.): 2.25