


Mounting Clamp for Medical-Grade Power Strips

Model IVSTAND-PSCC



Features:

- Mounts a medical-grade power strip or surge protector to an IV pole, crash cart or anywhere else needing portable power
- Patented antimicrobial coating is 99.9% effective in inhibiting the growth of bacteria
- Durable lightweight aluminum construction designed to stand up to heavy use
- Built-in drip shield protects outlets from liquid and other contaminants
- Mounting hardware included for quick, easy installation on IV pole, crash cart or bed
- Adjustable clamp fits poles ranging in diameter from .875 in. to 1.25 in.

Specifications:

- Color: White
- Antimicrobial: Yes
- Material of Construction: Metal
- Cable Management: Bar
- Shipping Dimensions (hwd / in.): 3.03 x 6.22 x 10.75
- Unit Weight (lbs.): 0.530