

Disposable 3-Ply Face Mask



Features:

- Adult Size (175mm x 95mm)
- Ear-loop Fastener
- Adjustable Nose Strip
- Available in BLUE or BLACK
- Cover the user's nose and mouth
- For Single or Disposable Use Only
- One Layer of Melt Blown Fabric encased in Two Layers of Non-Woven Fabric
- 50 Masks per Inner Box
- 40 Inner Box in a Master Carton
- Total 2,000 Masks/Carton
- Discounts available for orders of 12,000 Mask or more