



Medical Recliner/Transporter

Model RC301-A2P



Features:

- Hand-held remote
- Swing-away arms
- Fold-away side table (one)
- Push handle
- Retractable foot tray
- 4" twin wheel casters
- Wall-hugger mechanism (infinite recline positions)
- Quick-release seat
- All steel, powder-coated frame
- Patient weight capacity: 300 lbs.

Options:

- Articulating table
- Central locking caster system
- Trendelenburg
- Additional fold-away side table
- Arm cap
- Cupholder - On all tables or one of two tables
- Heat and massage
- Cal 133 fire barrier
- Vinyl upgrades:
 - Non standard vinyl*
 - Customer's own material*
 - PVC-free material*
- Dual-hand control
- Battery backup

Specifications:

- Back Height (from floor/seat) - 49.5"/28.75"
- Overall width (2 table down/up) - 32.75"/52"
- Seat Width - 22.5"
- Arm Height from Floor (front/back) - 29"/26"
- Seat Height/Depth - 24.5"/20.5"
- Upright Overall Depth (w/handle) - 40.5"
- Thermal Plastic Casters - 4"
- Weight Capacity - 300 lbs.