

"EASY" Medical Recliner/Transporter

RC301ATF | RC301PATF



Features:

- Infinite Recline (Trendelenburg)
- Fold-away side table
- Retractable foot tray
- 5" wheel casters
- Wall-hugger mechanism
- All steel, powder-coated frame
- Cal 133 fire barrier
- Patient weight capacity: 500 lbs.

Options:

- 2nd Side Table
- Cup Holder for Side Table
- Push Bar
- Manual or Powered Chairs



Specifications:

- 30" Width
- 44" Height (27" when fully reclined)
- 35" Depth (72" when fully reclined)
- 22" Seat Width
- 21" Seat Depth
- 27" Arm Height

