



Cast Iron IV Stands

NV-CIB-SSIV-FA



Features:

- Sturdy, Anti-Tipping, Cast Iron Base with a durable, black powder coat finish to allow for easy cleaning.
- 5-leg, low-profile base is heavy duty and offers extra stability.
- Heavy gauge steel pole for strength and durability.
- Easy-To-Use twist mechanism to raise and lower the pole.
- 3" Casters
- Aesthetically appealing look and overall design
- Well engineered product that will provide years of carefree service
- Height Range: Minimum 53.5" - Maximum 95" tall
- Base Diameter: Approximately 25"
- 5 Year Warranty
- Unit Weight: 37 lbs.



Options:

- Push Wheel
- Foley Hooks
- Steer Handle
- Power Strip