

Renaissance Privacy Screen Mobile Trolley

RPS-FG6



Features:

- Mounting on the Trolley chassis
- Max. 10-wing, shielding 300 cm
- 2 screen heights to choose from (145/165 cm)
- Trolley base 52 x 59 cm
- Individually breakable smooth running casters
- Easy to operate lifting mechanism
- Easy to disinfect/clean
- Movable when in "parked" position
- **Note:** Mounted on a trolley, the folding screen is mobile only. It is not possible to mount the device on a wall, device rail or column.



Specifications:

- Privacy Screen Height Options:
 - 145 cm (57.09 inches)
 - 165 cm (64.96 inches)
 - 185 cm (72.83 inches) *Not compatible with Trolley Mount*
- Trolley Mount: 52cm x 59cm (20.47inches x 23.23inches). Mobile only. Cannot be attached to wall or column.